


OxySafe[™] Options

(#1109-2522)



OxySafe should be placed as close to what you want to protect as possible.

Use OxySafe's special tubing to

Use OxySafe's special tubing to place OxySafe as close to the patient and as close to the oxygen supply as possible. We recommend the OxySafe Complete Patient Kit.

(# 1109-2522)

RECOMMENDED SET-UP

OxySafe Complete Patient Kit

Kit includes two OxySafes, one 2" tubing connector (#1109-2014), and one 2" nasal cannula (#1109-2023). Sold in bags of 25.

1109-2522

